

Confidence Down the Line

BD™ Select APS TSB

Improved DCM for your media fill applications:

- ✓ Manufactured in an Animal Origin-free, Antibiotic-free facility
- ✓ Mycoplasma Testing Upon Release
- ✓ Endotoxin Testing Upon Release



BD™ Media Fill Application Products

Select APS (Alternate Protein Source) TSB (Tryptic Soy Broth) Products

Cat. No.	Description	Packaging	Storage	Added Testing
214886	Select APS TSB Irradiated	10 Kg	2 to 25° C	Mycoplasma/Endotoxin
214887	Select APS TSB	10 Kg	2 to 25° C	Mycoplasma/Endotoxin
214889	Select APS TSB	500 g	2 to 25° C	Mycoplasma/Endotoxin

Standard TSB Products

Cat. No.	Description	Packaging	Storage
296264	Trypticase Soy Broth Sterile	500 g w/ outer bag	2 to 25° C
211823	Bacto Tryptic Soy Broth	10 Kg	2 to 25° C
211824	Tryptic Soy Broth	100g	2 to 25° C
211822	Tryptic Soy Broth	2 Kg	2 to 25° C
211768	Trypticase Soy Broth	500 g	2 to 25° C
257203	Tryptic Soy Broth	5 Liter Bag	2 to 8° C
257202	Tryptic Soy Broth	10 Liter Bag	2 to 8° C



Plated Media for Environmental Monitoring

Cat. No.	Description	Dish	Wrap	Packaging	Storage	SAL Claim
222205	Sterile Pack TSA 100mm 100 EA	100 mm	Double wrap	100 plates	2 to 8° C	10 ⁶
215343	IC-XT TSA 100MM 100 EA	100 mm	Triple wrap	100 plates	2 to 25° C	10 ⁶
215342	IC-XT TSA LP80 100MM 100 EA	100 mm	Triple wrap	100 plates	2 to 25° C	10 ⁶
292336	Sterile Pack Sabouraud Dextrose Agar	100 mm	Triple wrap	100 plates	2 to 8° C	10 ⁶
215276	Sterile Pack TSA LP80 Rodac Snap	RODAC™ Snap Lid	Double wrap	100 plates	2 to 8° C	10 ⁶
215285	Isolator Pack TSA LP80 Rodac Snap	RODAC™ Snap Lid	Triple wrap	100 plates	2 to 8° C	10 ⁶
215338	IC-XT TSA LP80 RODAC SNAP 100 EA	RODAC™ Snap Lid	Triple wrap	100 plates	2 to 25° C	10 ⁶
215341	IC-XT TSA LP80 RODAC 100 EA	RODAC™	Triple wrap	100 plates	2 to 25° C	10 ⁶
292337	Isolator Pack Sabouraud Dextrose Agar LP80	RODAC™	Triple wrap	100 plates	2 to 8° C	10 ⁶
215282	Sterile Pack Sabouraud Dextrose Agar LP80	RODAC™ Snap Lid	Double wrap	100 plates	2 to 8° C	10 ⁶

Bottled Media Products

Cat. No.	Bottle Description	Volume	Packaging	Storage
257543	TSB EtO Syrup	100 ml	20 ea double wrapped box	2 to 25° C
257544	FTM EtO Syrup	100 ml	21 ea double wrapped box	2 to 25° C
257097	FTM Special EtO Flip Cap	100 ml	44 ea single wrapped box tyvek	2 to 25° C
257307	FTM Special EtO Flip Cap	100 ml	44 ea single wrapped box tyvek	2 to 25° C
299107	TSB Serum	100 ml	10 ea	2 to 25° C
299108	FTM Serum	100 ml	10 ea	2 to 25° C
257485	FTM Wide Mouth	100 ml	10 ea	2 to 25° C
257486	TSB Wide Mouth	100 ml	10 ea	2 to 25° C
290652	Fluid A w/Septum	300 ml	10 ea	2 to 25° C
290662	Fluid D w/Septum	300 ml	10 ea	2 to 25° C
299111	Fluid A w/Septum	500 ml	10 ea	2 to 25° C
299110	Fluid D w/Septum	500 ml	10 ea	2 to 25° C



BD Life Sciences
 Diagnostic Systems
 7 Loveton Circle
 Sparks, MD 21152-0999 USA
 Tel: 800.638.8663
www.bd.com/ds